

The right pipette tip for your Transferpette®

Optimal results are achieved in combination with the BRAND micropipettes. Quality tips from BRAND provide a perfect fit for precise analysis. Pipette and filter tips from BRAND are tested for BRAND pipettes and for most pipette models from GILSON®, Thermo Fisher Scientific FINNPIPETTE®, Eppendorf®, and Sartorius®.

The 5 ml tip can be used with the pipette models from models from BRAND and Thermo Fisher Scientific FINNPIPETTE®. The 1250 µl and 10 ml tips are suitable for BRAND and Eppendorf®.

Note: Pipette shafts may be subject to modifications. Compatibility depends, among other things, on the manufacturer, pipette type, serial number, and manufacturing date. We recommend using the free sample bag to check tip compatibility.



Transferpette® S

Transferpette® electronic

Pipette tips Filter tips Volume range	Transferpette® single-channel Nominal volume													Transferpette® multi-channel Nominal volume					
	1 µl	2,5 µl	10 µl	20 µl	25 µl	50 µl	100 µl	200 µl	500 µl	1000 µl	2 ml	5 ml	10 ml	10 µl	50 µl	100 µl	200 µl	300 µl	1250 µl
0.1 – 20 µl	✓	✓ ₂₎	✓ ₂₎											✓					
0.5 – 20 µl		✓	✓											✓					
1 – 50 µl		✓ ₂₎	✓											✓					
2 – 200 µl ¹⁾				✓	✓	✓	✓	✓							✓	✓	✓	✓	
5 – 300 µl				✓	✓	✓	✓	✓							✓	✓	✓	✓	
50 – 1000 µl ¹⁾									✓	✓									
50 – 1250 µl ³⁾										✓ ₂₎									✓
50 – 1250 µl XL ³⁾										✓ ₂₎									✓ ₂₎
0.5 – 5 ml											✓	✓							
1 – 10 ml													✓						
0.1 – 1 µl	✓	✓ ₂₎	✓ ₂₎											✓					
0.5 – 10 µl		✓	✓											✓					
1 – 20 µl			✓											✓					
2 – 20 µl				✓	✓	✓	✓	✓							✓	✓	✓		
5 – 100 µl				✓	✓	✓	✓	✓							✓	✓	✓	✓	
5 – 200 µl								✓							✓	✓	✓	✓	
50 – 1000 µl									✓	✓									
50 – 1250 µl ³⁾										✓ ₂₎									✓
50 – 1250 µl XL ³⁾										✓ ₂₎									✓ ₂₎

✓ = Tip volume less than pipette's nominal volume

¹⁾ Bulk tips are yellow or blue colored; racked tips are colorless in a yellow-colored or blue-colored tip tray

²⁾ Calibration and possible adjustment needed

³⁾ 1250 µl = Length 86 mm, 1250 µl XL = Length 102 mm

